






TOPAIR
CLEAN AIR SOLUTIONS

IPMS - Integrated Particle Monitoring System

The IPMS measures particle concentration in real time, monitoring the cleanliness level 24/7. The clear display indicates whether the cleanliness level complies with ISO 5 and alerts when the workspace has been contaminated or needs to be serviced.


The IPMS smart algorithm processes data provided by an advanced built-in sensor and displays a reliable outcome. IPMS provides critical real-time information on the air quality within the workspace, prevents cross contamination, and ensures constant compliance with relevant standards.

The IPMS is an unprecedented solution that can transform the future of biological processes in the industry.



GREEN LIGHT

Chamber environment is equal or better than the IOS-5 requirement within the workspace



RED LIGHT

Chamber environment is lower than the IOS-5 requirement within the workspace