



TOPAIR
CLEAR AIR SOLUTIONS

Opioid Evidence Handling Workstation Kit

TopAir’s Opioid Evidence Handling Workstation offer dual functionality by providing both fume containment and precise weighing capabilities.

The stable cabinets are designed for accurate measurements of powders and chemicals by delivering a weight precision up to four decimal places.

- Dual ULPA filter The ULPA filter includes a main filter and a secondary back up filter. It removes powders and particles with an efficiency of 99.9995 @ 0.3µm and provides reassurance to the to the lab staff.
- Vacuum system for efficient residue collection The vacuum system is an excellent instrument to collect residues of powders from the working surface and channel them to the ULPA filter.
- Filter integrity test after production (D.O.P) This test provides reassurance to lab staff about air quality and contamination control in the laboratory environment

Spec/ Model	CF-060-OP	CF-090-OP	CF-120-OP	CF-150-OP	CF-180-OP
Opioid kit for Duct-less Fume Hood	Opioid kit for Duct-less Fume Hood	Opioid kit for Duct-less Fume Hood	Opioid kit for Duct-less Fume Hood	Opioid kit for Duct-less Fume Hood	Opioid kit for Duct-less Fume Hood

