

Weighing Station Kit

Part Number : CF-W

TopAir's Ductless Fume Hoods offer dual functionality by providing both fume containment and precise weighing capabilities.

The stable cabinets are designed for accurate measurements of powders and chemicals by delivering a weight precision up to four decimal places.

For enhanced weighing performance, an optional "weighing station kit" is available. This kit includes:

- HEPA filter
- Vacuum system for efficient residue collection

The HEPA filter removes powders and particles with an efficiency of 99.9995 @ 0.3 μ m and provides reassurance to the lab staff. The vacuum system is an excellent instrument to collect residues of powders from the working surface and channel them to the HEPA filter.

Structure:

- Unit is based on the standard ductless unit (see pages 2-3)



Clear Back for Ductless Fume Hood (Educational)

Part number: CF-CB

TopAir's Educational Ductless Fume Hood with 360° glass provides an effective tool for educational purposes and a safe environment for lab staff working with acids and harsh chemicals.

Structure:

- The unit is based on the standard ductless unit (see pages 2-3)
- The unit includes a clear back panel made out of polycarbonate.

