

# Ductless Fume Hood Versions

## XL Double Filter for Ductless Fume Hood

**Part number: CF-XL**

The Ductless Fume Hood XL unit builds upon the original unit (pages 2-3) by offering increased versatility with space for two filters. This allows you to customize filtration based on your specific needs. Filters are sold separately.

The XL model is compatible with a wide range of filters, including HEPA, ULPA, carbon filters for solvents, acids, bases, and formaldehydes, as well as multi-gas filters.

Our team is happy to assist you with selecting the most appropriate filters to ensure optimal protection for your application

Model	Height
CF-XL	1,450 mm (instead of 1,223 of the original unit) 57.1" (instead of 48")

